

The gold standard of removable dental appliances care - Dr. Mark's HyGenie® ergonomic oral hygiene device

This spring, while attending the 38th International Dental Exhibition (IDS), Cologne, Germany, I came across Dr. Mark's HyGenie® ergonomic oral hygiene device as a first innovation. Despite innovations dedicated to mobile dental devices, their care has remained unchanged.

Dr. Mark's HyGenie®, a simple solution to a common problem, has been created by an Australian, Dr. Mark Woakpoon, BDS, who, after four years and thousands of hours of work dedicated to the creation process, has made the device available to the consumers as part of a project supported by the Australian Government Department of Industry, Innovation and Science based on marketing grant under the Entrepreneurs' Programme. In addition to Dr. Mark's HyGenie® device, Dr. Woakpoon has also developed Dr. Mark's SureGrip SureGrip® and Dr. Mark's DentalFresh.

Dr. Mark's HyGenie device is multifunctional, serves all types of removable dental appliances, cleans, stores and protects them.

The device made of durable materials lasts up to two years in daily use, proving its efficiency, accessibility and elegance. Each component was designed, and each material was chosen for the purpose of high-performance care of all types of removable dental devices.

The ergonomic device fits all ages, has a set of *clustered brushes* that rotate effortlessly in opposite directions. After inserting the removable dental device into the oral hygiene device, apply two to three drops of DentalFresh, pass the device under a stream of running water, make about 15-20 turns, and then flush it.

The rotating action causes the brushes to clean all surfaces of the removable dental device simultaneously.

The cleaned dental device can be worn immediately or stored in the oral hygiene device.

Dr. Mark's DentalFresh is a product designed to keep your dental device clean and fresh. The solution is a formula with pure peppermint oil, paraben- and alcohol-free, designed for daily use with a pH set to 5.0, which dissolves calculus accumulation, without affecting the dental apparatus.

I use and recommend Dr. Mark's HyGenie® multifunctional ergonomic device for clean, store, and protect purposes, to all my patients who carry different removable dental devices from sports mouthguards, orthodontic aligners, night guards, splints, retainers, to full dentures, overdentures, and partial dentures.


Let me also take this opportunity to warmly recommend it to you.



Dr. Mark's HyGenie device (Dr. Mark's HyGenie Pty Ltd., Perth, Wa 6000, Australia) <https://drmarkshygenie.com/>.



Dr. Mark's DentalFresh (Dr. Mark's HyGenie Pty Ltd., Perth, Wa 6000, Australia) <https://drmarkshygenie.com/>.



Dr. Mark's SureGrip Denture Brush (Dr. Mark's HyGenie Pty Ltd., Perth, Wa 6000, Australia) <https://drmarkshygenie.com/>.

Florin - Eugen Constantinescu
DMD, PhD Student
Editorial Director, Product News

DOI: [https://doi.org/10.25241/stomaeduj.2019.6\(3\).prodnews.1](https://doi.org/10.25241/stomaeduj.2019.6(3).prodnews.1)